


Specification Overview



Instrument type	Benchtop
Dimensions	Approx. H397 × W168 × D488 mm
Weight	Approx. 22 kg
Detection System	Chemiluminescence or Fluorescence

Layout (Reagent Cartridge)



All Process In Tip technology

Multiplex
Quantitative
Highly sensitive
Fully Automated

Micro Plate → 1 mm beads

LuBEA System

Produced by Precision System Science Co., Ltd.

- 20 Beads in Tip
 - Liner Prefilled Reagent Cartridge
 - 1 mm immuno assay
 - 1 mm DNA Typing
 - Chemiluminescence
 - Fluorescence
-) 2 Types Detection method



● For customers in Asia / Pacific
Precision System Science Co., Ltd.
 88 Kamihongou, Matsudo-shi, Chiba 271-0064, Japan
 Tel: +81-47-303-4801 Fax: +81-47-303-4811
 URL: <http://www.pss.co.jp>
 E-mail: service@pss.co.jp

● For customers in North / South America
Precision System Science USA, Inc.
 Tel: +1 (925) 960-9180 / FAX: +1 (925) 960-9184
 E-mail: contact@pssbio.com

● For customers in Europe / Africa / Middle East
Precision System Science Europe GmbH
 Tel: +49 (0) 6131 6966 468 / FAX: +49 (0) 6131 6966 469
 E-mail: contact-pss@pss.co.jp

○ The performance, specifications and appearance of products described in this catalogue are subject to change without prior notice.
 ○ The information in this catalogue is current as of January, 2015.

Device Overview

LuBEA is a small, automated detection device which uses the multi-item detection tool "BIST" (Beads array In Single Tip).

An original PSS technology, BIST allows the building of an automatic reaction protocol, and when further combined with pre-filled reagent, enables automation of reaction, washing, and detection.

With the combination of BIST, test sample, and pre-filled reagent, detection can be made without any manual operation.

The automated device is easy to use for a variety of fields, including polymorphic DNA-based tailor made direction, clinical testing, and food.

※ We are currently developing a detection unit that combines both chemiluminescence and fluorescence to allow selection of multiple items to be simultaneously detected using one BIST.


About BIST

BIST is an analytical signal detection tool loaded with 1mm-diameter beads which each have different functions. Beads with differing DNA fragments and antibodies are lined up in the capillary, and each functions the same as an individual reaction chamber. The tool allows a wide array of tests through the insertion of various bead types.

By shaping the capillary as a dispensing tip and combining it with the LuBEA automated system, we enable simultaneous measurement of multiple items and complete automation of the reaction process.

When compared with previous methods requiring

separate detection for each reaction, BIST allows for a smaller sample size and reagent volume needed to obtain results. The reaction process is also sped up through greater liquid mixing efficiency.




Test items requiring Multiplex (Ex.)

Body testing for health guidance (obesity factors), talent analysis, other applications	β 3AR, UCP1, AGT235, MTHFR	Whole blood, serum, plasma
Allergic reaction tests	Cedar, cypress, hogweed, house dust	
Cancer markers	AFP (Liver cancer), PSA (prostate cancer), P53 antibody	
Cancer-related tests	HLA Typing, sugar chain detection	
Cardiac disease marker	Myoglobin, troponin, PCT	
Clinical testing items	TSH, Free T3, Free T4, AFP, albumin	
Autoimmune diseases		
Pharmacogenomics	IL28, CYP2C19, EGFR (Iressa target agent)	
Immunodeficiency disease	CD3 (T lymphocyte), CD19 (B lymphocyte), CD45 (leukocyte)	
Infant mass screening	TSH, 17-OHP, Free T4	
Genetically modified food (GMO) proteins	Herbicide-resistant enzyme PAT, CP4-EPSPS	Foods (corn, soybeans)
Genetically modified food (GMO) genes	Specific gene detection	


BIST for ELISA Immunoassay

This is a BIST for immunoassay that allows ELISA by lining up beads affixed with antibodies and antigens in the capillary. Results can be gained with certainty even for low concentration samples (such as antibodies and antigens) by efficiently mixing reaction fluid and establishing a high sensitivity detection system.



BIST for DNA Typing

DNA testing requires numerous steps, including nucleic acid extraction from the sample, PCR, SNP detection reaction, and BIST reaction. We have developed an instrument that performs all processes from nucleic acid extraction to BIST detection for a stable testing process that avoids the complexity of manual operation.



Food allergens	Wheat, buckwheat, peanuts, eggs, milk, shrimp, crab	Food
Crop production location identification, variety identification	Tuna variety identification, cashmere identification	Meat, hair, etc.
Food poisoning	Salmonella (bacteria)	Stool, food
	Staphylococcus aureus (bacteria, toxins: enterotoxins)	
	Vibrio parahaemolyticus (bacteria)	
	Bacillus cereus (bacteria, toxins)	
Legal communicable diseases	Cholera, dysentery, typhoid	Water
	Diarrheal viruses	
Enterohemorrhagic E. coli	O157, O26	
Communicable diseases	Mumps, EV virus	
Communicable diseases	Legionella	